

### KINGSTAR



www.antigentestdevice.com



# Type IIR Disposable Surgical Face Mask

Find a huge selection of Type IIR
Disposable Surgical Face Mask from
China at KINGSTAR INC. For clinical
medical personnel to wear during
operation. Type IIR Disposable Surgical
Face Mask covers the user's mouth,
nose, and jaw, providing a physical
barrier to prevent the direct passage of
pathogens,microorganisms, body fluids,
and particles.

KINGSTAR INC is Type IIR Disposable Surgical Face Mask manufacturers and suppliers in China who can wholesale Type IIR Disposable Surgical Face Mask.

### **Product Introduction**

Type IIR Disposable Surgical face mask is a type of medical face mask which has higher protection than normal medical face mask. This type of face mask is more fit for nowadays environment to protect yourselves.

# **Specification**

Packing Specifications: 10pcs/pkg, 5pkgs/box, 48 bxs./CTN

Box Dimensions: 19.5×10×8.5cm Carton Dimensions: 41.5×40.5×53cm

Weight: 11.5 kg











# **Product Feature and Application**

Type IIR Disposable Surgical face mask is flat ear loop plane-type with size 17.5cm×9.5cm. Bacterial filtration efficiency is not less than 98%



# Type IIR Disposable Surgical face mask











3 PLY

Breathable

Protection

Disposable

### **Product Details**

Type IIR Disposable Surgical face mask has 3 layers.

First layer——Polypropylene spunbonded non-woven fabric

Second layer——Polypropylene melt-blown non-woven fabric

Third layer——Polypropylene spunbonded non-woven fabric


Nose clip——Polyethylene and iron wire

Rubber bands——Nylon / Polyurethane composites



# Three-layer Disposable medical face mask

B.F.E. (Bacterial filtration efficiency) ≥ 98%




# Built-in nose clip



Three-dimensional nose bridge clip for easy fit



# Extra-soft earloops



# Comfortable to wear no hurt to ears

# Warning and precaution

- 1. Please do not use the mask when the mask is contaminated or the whole mask is damaged and deformed.
- 2. This mask is disposable only and cannot be reused.
- 3. The side with the nose clip is facing outward when wearing.

# Storage and stability


Type IIR Disposable Surgical face mask is recommended to store in a clean, dry, ventilated and non-corrosive gas environment, away from ignition sources or combustible materials.

The following is the certificates of Type IIR Disposable Surgical face mask.











# INGEER CERTIFICATION ASSESSMENT SERVICES



# MANAGEMENT SYSTEM CERTIFICATE

Address Middle of Naming Road, Langue (Natrict, Chuale

NED-ver East of Naming Stead, 198-of of Benglin Road, Smith of Annying Soud, Namits of Singling Soud, Langue (Nativa), Churchen City, Anbuij Paneirace, China

Certificate No.: 117 20 QOM 0146 ROM
This is to certify the quality management systems of

Unified Social Credit Code 91341102MA2RRK2RON

has been assessed and registered as meeting the requirements of ISO13485:2016 <Medical device-Quality management systems - Requirements for regulatory purposes>

City, Aphal Province, China

Production of Disposable Medical Mask and Disposable Surgical Mask

Signed by:



First Certification: 11 Sep. 2020 Expiry Date: 10 Sep. 2023

# Shanghai Ingeer Certification Assessment Co.,Ltd.

Certification and Accreditation Administration of PRC:CNCA-R-2003-117
Tel: 400-182-9001/+86 21-51114700

Web: www.icas.org.cn

Add: Room 801;HuaDing Mansion, 2368# West Zhongshan Rd., Xuhui District, Shanghai, China, 200235



tr Sunsilance Audit Second Sunsilance Audit

The ownership of the certificate belongs to Shanghai Ingeer Certification Assessment Co., Ltd. The information & validation of this certificate can be checked on the CNCA website. WOW. CNCA.GOV. CN and ICAS website, or by calling ICAS's clients services Dept. The certificate is only valid when used together with related permits when appropriate. If the organization can't effectively maintain the above management system, ICAS has the right to withdraw the qualification certificate.